


Thermaloy® (Thermal Activated Nickel-Titanium) Preformed Natural and Ideal Arches


- Pliable at room temperature and superelastic at body temperature
- Same elastic working range as superelastic with lower forces
- Can be used earlier in treatment to provide tongue control
- Can be used as an initial arch wire for treatment
- Easier to ligate in severely maligned cases and impacted canines
- Increases patient comfort
- Permanent midline etchmarks: Three lines for a maxillary arch, one line for a mandibular arch

PREFORMED NATURAL ARCHES



Inch	Size (mm)	Maxillary	Mandibular	Package Contains
ROUND				
.013	(0.330)	A07310	A07330	10
.014	(0.356)	A07300	A07320	10
.015	(0.381)	A07311	A07331	10
.016	(0.406)	A07301	A07321	10
.017	(0.432)	A07312	A07332	10
.018	(0.457)	A07302	A07322	10
SQUARE				
.016 x .016	(0.406 x .0406)	A07303	A07323	10
RECTANGULAR				
.016 x .022	(0.406 x 0.559)	A07304	A07324	10
.017 x .025	(0.432 x 0.635)	A07305	A07325	10
.018 x .025	(0.457 x 0.635)	A07306	A07326	10
.019 x .025	(0.483 x 0.635)	A07307	A07327	10
.021 x .025	(0.533 x 0.635)	A07308	A07328	10

PREFORMED IDEAL ARCHES



Inch	Size (mm)	Maxillary	Mandibular	Package Contains
ROUND				
.013	(0.330)	A07360	A07380	10
.014	(0.356)	A07350	A07370	10
.015	(0.381)	A07361	A07381	10
.016	(0.406)	A07351	A07371	10
.017	(0.432)	A07362	A07382	10
.018	(0.457)	A07352	A07372	10
SQUARE				
.016 x .016	(0.406 x .0406)	A07353	A07373	10
RECTANGULAR				
.016 x .022	(0.406 x 0.559)	A07354	A07374	10
.017 x .025	(0.432 x 0.635)	A07355	A07375	10
.018 x .025	(0.457 x 0.635)	A07356	A07376	10
.019 x .025	(0.483 x 0.635)	A07357	A07377	10
.021 x .025	(0.533 x 0.635)	A07358	A07378	10